

POST OP FOOTWEAR

OAPL SQUARE TOE SHOE



OAPL POST OP SHOE



OAPL BUNION BOOT



OAPL ADL SHOE



OAPL CAST SHOE



OAPL FOREFOOT UNLOADER



OAPL SQUARE TOE SHOE*

A multipurpose post op shoe designed to provide protection and accommodate bulky dressings or bandages. The square toe acts as a bumper around the toes and easy to adjust Velcro tabs keeps the shoe secure.

| StockCode | Size | Length |
|-----------|------------|--------|
| SQTP | Paediatric | 205mm |
| SQTXXS | XX/Small | 230mm |
| SQTXS | X/Small | 245mm |
| SQTS | Small | 260mm |
| SQTM | Medium | 275mm |
| SQTL | Large | 290mm |
| SQTXL | X/Large | 310mm |

OAPL POST OP SHOE

Designed to wear after soft tissue procedures or post op care. Adjustable straps accommodate the bulkiest dressings. Non-skid sole absorbs shock during ambulating.

| StockCode | Size | Length |
|-----------|----------|--------|
| 11082C | XX/Small | 240mm |
| 11082D | X/Small | 255mm |
| 11082F | Small | 265mm |
| 11082G | Medium | 285mm |
| 11082H | Large | 300mm |
| 11082I | X/Large | 310mm |

OAPL BUNION BOOT

Wear post operatively after surgeries such as bunionectomies and corrective osteotomies. Also suitable for un loading the toes and forefoot plantar wounds.

| StockCode | Size | Length |
|-----------|--------|--------|
| 11082BUN | Small | 165mm |
| 11083BUN | Medium | 175mm |
| 11084BUN | Large | 185mm |

OAPL ADL SHOE

Suitable for treatment of the high risk and diabetic foot. The ADL boot comes with an OAPL HEX OFFLOADER inlay to allow for decreasing pressure through wounds and calloused areas. Wear with additional forefoot cap for added protection. Forefoot cap sold separately.

| StockCode | Size | Length |
|-----------|---------|--------|
| ADLS | Small | 240mm |
| ADLM | Medium | 260mm |
| ADLL | Large | 285mm |
| ADLXL | X/Large | 305mm |

OAPL CAST SHOE

A canvas sandal with a rubber rocker sole and Velcro closure. To protect plaster and fibreglass casts from dirt and moisture. Provides secure non-slip footing.



| StockCode | Size | Length |
|-----------|------------|--------|
| 11079A | Paediatric | 180mm |
| 11079 | Small | 240mm |
| 11080 | Medium | 280mm |
| 11081 | Large | 320mm |
| 11082 | X/Large | 360mm |

OAPL FOREFOOT UNLOADER

Suitable for unloading metatarsal heads, the toes and the forefoot. For use with diabetic feet and post operatively. Very adjustable across the forefoot and dorsum of the foot for bulky dressings. Breathable mesh material for added comfort. Tension firmly for optimal stability comes with optional comfort inlay.

| StockCode | Size | Length |
|-----------|---------|--------|
| FFUXS | X/Small | 240mm |
| FFUS | Small | 255mm |
| FFUM | Medium | 270mm |
| FFUL | Large | 285mm |
| FFUXL | X/Large | 300mm |

OAPL HEX OFFLOADER INSERT



10mm tri layer foam rubber insert with hexagon shaped pull out pieces, remove and provide offloading to the diabetic or sensitive foot.

| StockCode | Size | Length |
|-----------|---------|--------|
| SQTHEXS | Small | 260mm |
| SQTHEXM | Medium | 275mm |
| SQTHEXL | Large | 290mm |
| SQTHEXXL | X/Large | 310mm |

*Shoe is compatible with the OAPL HEX OFFLOADER INSERT. Prices do not include GST.

Orthopaedic Appliances Pty Ltd reserves the right to make any amendments and/or variations to this product or packaging without notice.